

Skin Tightening

How does Skin Tightening work?

ST stands for Skin Tightening and it uses near-infrared pulsed light technology to treat wrinkles and skin laxity through deep dermal heating that selectively targets collagen, protein and fat without damage to the epidermis. It heats the inner layer or dermis of the skin, inducing the formation of collagen, essential for a firm, young look. It affects skin tightening, toning and texturing.

Skin Tightening treatment is different from any other aesthetic device for collagen stimulation because the 780-1000nm wavelengths have been adjusted to target the collagen that they are trying to rebuild, without it being painful. Other skin tightening machines have a lot of pain associated with them because they haven't been able to filter out the high water adsorption wavelengths that cause pain.

How soon can results be seen?

Results can be seen immediately and will steadily improve over the following weeks and months.

Is there any discoloration after treatment?

Not generally. The advantage of pulsed light assures safe, effective treatment by targeting the dermis. This results in little or no skin damage or discomfort.

Is there any downtime with this treatment?

No. Clients experience little or no pain with the treatment, and can return to their normal activities immediately.